

SILICONE NURSE WATCH - round - purple



Product Code	20403
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	Not required
CH-REP	Not required
EAN/UPC	8023279204032

Description

Nurse Watch, round

Colour: purple

Silicone nurse watches, provided with a safety pin.
Suitable to measure patient pulse, beat and sterilization procedure timing.
Conforming to procedures not allowing nurses to wear wrist watches.

Watch case 35 mm, watch length 85 mm.

Shelf life of battery: 12-15 months from first use.

Manual: GB, FR, IT, ES, PT, DE, GR, PL , Arabic.