

SILICONE NURSE WATCH - square - hearts



Product Code	20461
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	Not required
CH-REP	Not required
EAN/UPC	8023279204612

Description

Nurse Watch, squared

Colour: hearts fantasy

Silicone nurse watches, provided with a safety pin.
Suitable to measure patient pulse, beat and sterilization procedure timing.
Conforming to procedures not allowing nurses to wear wrist watches.

Watch case 35 mm, watch length 85 mm.

Shelf life of battery: 12-15 months from first use.

Manual: GB, FR, IT, ES, PT, DE, GR, PL, Arabic.