



| | |
|---------------|-------------------|
| Product Code | 20464 |
| Unit of sale | 1 pc |
| Minimum order | 1 |
| Type | No medical device |
| Class | |

| | |
|---------|---------------|
| UK-REP | Not required |
| CH-REP | Not required |
| EAN/UPC | 8023279204643 |

Description

Nurse Watch, squared

Colour: roses fantasy

Silicone nurse watches, provided with a safety pin.
Suitable to measure patient pulse, beat and sterilization procedure timing.
Conforming to procedures not allowing nurses to wear wrist watches.

Watch case 35 mm, watch length 85 mm.

Shelf life of battery: 12-15 months from first use.

Manual: GB, FR, IT, ES, PT, DE, GR, PL, Arabic.