

S/S BANDAGE SCISSORS - purple ring - 14 cm



Product Code	20502
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
CH-REP	Yes

RDM (NSIS)	2775702
CND	L01040199
EAN/UPC	8023279205022
GTIN	08023279205022
GMDN	13481

Description

Stainless steel nurse bandage scissor - 14 cm

- Model: Lilac ring

Powder coated instruments are not autoclavable, can be sterilized in water at maximum 110°C