

# S/S STRAIGHT SCISSORS with belt clip - blunt/sharp - 14 cm



|               |                |
|---------------|----------------|
| Product Code  | 20520          |
| Unit of sale  | 1 pc           |
| Minimum order | 1              |
| Type          | Medical device |
| Class         | I              |

|        |     |
|--------|-----|
| UK-REP | Yes |
| CH-REP | Yes |

|            |                |
|------------|----------------|
| RDM (NSIS) | 2775713        |
| CND        | L01040199      |
| EAN/UPC    | 8023279205206  |
| GTIN       | 08023279205206 |
| GMDN       | 13481          |

## Description

Stainless steel nurse straight scissor - 14 cm

- Model: Standard with belt clip

Powder coated instruments are not autoclavable, can be sterilized in water at maximum 110°C