

S/S STRAIGHT SCISSORS - leopard fantasy - blunt/sharp - 14 cm



Product Code	20531
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
CH-REP	Yes

RDM (NSIS)	2775713
CND	L01040199
EAN/UPC	8023279205312
GTIN	08023279205312
GMDN	13481

Description

Stainless steel nurse straight scissor - 14 cm

• Model: Leopard

Powder coated instruments are not autoclavable, can be sterilized in water at maximum 110°C