

# GISAFE 98% FILTERING SURGEON MASK 3 PLY type IIR with loops - pediatric - light blue - flowpack



EN 14683:2019 + AC:2019  
BFE > 98%

Product Code	20680
Unit of sale	box of 10 pcs
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	No
RDM (NSIS)	2512790
CND	T020604
EAN/UPC	4897121890100
PARAF	981067630

## Description

98% filtering 3 ply masks - TYPE IIR with earloops

- Pediatric masks
- Flowpack of 10
- Coloured: light blue

Lightweight and comfortable surgical masks suitable for hospital, dental and outpatient use.

- conformable nose band for anatomical fitting
- fiber glass free, filtration efficiency BFE =98%
- non woven material, latex free
- 1st ply, outside 28 gsm non woven
- 2nd ply, 25 gsm melt blown filter
- 3rd ply, inside 25 gsm non woven

Multilingual box: GB, FR, IT, ES, PT, DE, SE, GR, RO, Arabic.

## Technical Specifications

BFE% = 98  
Differential pressure (Pa/cm<sup>2</sup>) < 60  
CfU/g < 30  
Size cm: 17.5 x 9.5 (adult)  
EN 14683:2019 + AC:2019 conform; MDR 2017/745 Class I