



Product Code	20716
Unit of sale	box of 20 pcs
Minimum order	1
Type	PPE
Class	III

UK-REP	Not required
EAN/UPC	4046719563821
PARAF	987476102

Description

3M™ AURA™ FFP3 RESPIRATOR 1863+ IIR

- Filtering efficiency: FFP3 NRD >99%

The 3M™ Aura™ Particulate Respirators provide effective respiratory protection for use in medical environments where health care workers are exposed to airborne dust particles, non-volatile liquid particles and bio-aerosols. Limiting the transmission of infective agents from staff to patients and suitable for use during surgical procedures and certain other medical procedures.

Resistance to penetration of liquid splashes.

Key Features

- tested to EN 14683:2005 (Surgical masks - Requirements and test methods) and EN 149:2001+A1:2009 (Respiratory protective devices - Filtering half masks to protect against particles - Requirements, testing, marking)
- CE mark according to PPE Regulation 2016/425 and in compliance to EN 14683:2005 Type IIR as surgical mask according to Medical Device Directive 93/42/CEE
- foldable, proprietary 3-panel design allows for greater facial movement and wearer comfort and easy storage when not in use
- low breathing resistance filter technology gives effective filtration with low breathing resistance for consistent high quality performance
- sculpted nose panel helps conform to the nose and contours of the face and helps to improve compatibility with 3M eyewear
- innovative chin tab for easy positioning, donning and adjustment
- individual hygienic packaging protects the respirator from contamination before use
- large, soft nosefoam is comfortable on the skin
- headbands accommodate the shape of your neck, face and head
- outer cover provides resistance to fluid splashes
- coloured headbands for easy identification: blue for FFP2 and red for FFP3

Multilingual box: GB, FR, IT, ES, PT, DE, NL, SE, DK, NO, FI, IS, GR, PL, HU, CZ, SK, SI, EE, LV, LT, RO, RU, UA, BG, TR.