

# MYCROCLEAN KID REUSABLE SURGICAL MASK - BFE 99.8% - blue-blue



Product Code	20776
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	No
RDM (NSIS)	1964086
CND	T0399
EAN/UPC	8023279207767
PARAF	981497391

## Description

Reusable surgical protective masks classified Type II according to the UNI EN 14683:2019 standard.  
Made of polyester 99% and carbon fiber fabric 1%.  
Materials are of EEC preferential origin. Patented design.  
Washable, by hand or in the washing machine at 60°C/140°F.

2020V3 masks for health professional and outdoor/indoor use  
One layer. Washable up to 500 times at 60°C/140°F.  
Possibility to customize with 1 colour logo (minimum 1.000 pcs)

- Size: Kid (age 6-9)
- Colour: blue
- Hem: blue
- Nose clip: NO

### PFE & BFE values

Particle Filtration (%) 0.3 µm: 90% - 0.5 µm: 95%.

Certified non-cytotoxic by tests carried out according to UNI EN ISO 10993-5. BFE = 99.8% (from laboratory tests carried out in accordance with UNI EN 14683:2019).

### Biological load

UFC < 30 UFC (from laboratory tests carried out in accordance with UNI EN 14683:2019). Air permeability (cm<sup>2</sup>/cm<sup>3</sup>/sec) ASTM E96 (also known as ISO-11092 or Ret or Hohenstein test) = 7,5.

All masks have multilingual instructions: GB, FR, IT, ES, PT, DE, GR, SE, PL, RO, Arabic.

Made in Italy.