



Product Code	20833
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
CH-REP	Yes

RDM (NSIS)	2837034
CND	T0207
EAN/UPC	8023279208337
PARAF	980297030
GTIN	08023279208337
GMDN	38260

Description

Wide range of patterns to meet all healthcare professionals needs.
Light and comfortable, provided with adjustable laces.

Available in two sizes: medium size (M), suggested for men and women with short or medium-long hair; large size (L), suggested for women with long, curly hair.

100% high quality pure cotton, highly breathable surgical caps. MD Class I.