

# FOLEY 2-WAY SILICONE CATHETER ch/fr 18 - ball 30 ml - sterile



Product Code	22097
Unit of sale	box of 10 pcs
Minimum order	1
Type	Medical device
Class	II A

UK-REP	No
CH-REP	Yes
RDM (NSIS)	2842241
CND	U010201
EAN/UPC	8023279220971
PARAF	977809944

## Description

100% SILICONE STERILE FOLEY CATHETERS - SINGLE USE - DEHP FREE  
Two-way catheters - balloon 30 ml - length 40 cm - hard valve

- ch/fr 18

Multilingual box and manual: GB, FR, IT, ES, PT, DE, NL, PL, CZ, SK, HU, EE, LT, LV, DK, SE, FI, RO, BG, GR, HR, SI