

# SURGICAL GOWNS 40 g/m2 140x160 cm - size XXL - non sterile



Product Code	21043
Unit of sale	box of 100 pcs
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	No

RDM (NSIS)	2613065
CND	T020401
EAN/UPC	6973866552034
UMDNS	11901

## Description

Surgical gowns 40g - non sterile

- Size: XXL
- Length: 140 cm
- Width: 160 cm
- Length of sleeve: 61 cm
- Packing size: 66 x 35 x 43 cm
- Box of 100 (10 pouches of 10 pieces)

Disposable surgical gowns that guarantee comfort and protection.

- fluid repellent
- abrasion resistant
- ultrasonic sewing
- book folding
- knitted cuff/velcro at neck
- anti microorganism and antibody fluids
- antistatic treatment
- great air permeability
- high strength and reliability
- set-in sleeve

Material: SMMS (spunbond + 2 meltblown + spunbond non wovens).

Weight per piece: 120 ± 2g

Multilingual instructions: GB, FR, IT, ES, PT, DE, GR, SE, RO, PL, Arabic.