


REINFORCED SURGICAL GOWNS 50 g/m2 120x140 cm

- size M - sterile



| | |
|---------------|----------------|
| Product Code | 21052 |
| Unit of sale | box of 50 pcs |
| Minimum order | 1 |
| Type | Medical device |
| Class | I S |

| | |
|--------|----|
| UK-REP | No |
| CH-REP | No |

| | |
|------------|---------------|
| RDM (NSIS) | 2728032 |
| CND | T020402 |
| EAN/UPC | 6973866550160 |
| UMDNS | 11897 |
| PARAF | 983837600 |

Description

Reinforced surgical gowns 50g - sterile

- Size: M
- Length: 120 cm
- Width: 140 cm
- Length of sleeve: 58 cm
- Packing size: 66 x 37 x 40 cm

Disposable surgical gowns that guarantee comfort and protection.

- fluid repellent
- abrasion resistant
- ultrasonic sewing
- book folding
- knitted cuff/velcro at neck
- anti microorganism and antibody fluids
- antistatic treatment
- great air permeability
- high strength and reliability
- reinforced pad 38 g (PP + PE) on the chest (50x70 cm) and sleeve (30 cm)
- set-in sleeve

Material: SMMS (spunbond + 2 meltblown + spunbond non wovens).

Weight per piece: 120 ± 2g

Sterile gowns are individually packed.

Sterilized by oxide Ethylene.

Multilingual instructions: GB, FR, IT, ES, PT, DE, GR, SE, RO, PL, Arabic.