


REINFORCED SURGICAL GOWNS 50 g/m2 120x150 cm - size L - sterile



Product Code	21053
Unit of sale	box of 50 pcs
Minimum order	1
Type	Medical device
Class	I S

UK-REP	No
CH-REP	No

RDM (NSIS)	2728032
CND	T020402
EAN/UPC	6973866550177
UMDNS	11897
PARAF	983837612

Description

Reinforced surgical gowns 50g - sterile

- Size: L
- Length: 120 cm
- Width: 150 cm
- Length of sleeve: 59 cm
- Packing size: 66 x 37 x 40 cm

Disposable surgical gowns that guarantee comfort and protection.

- fluid repellent
- abrasion resistant
- ultrasonic sewing
- book folding
- knitted cuff/velcro at neck
- anti microorganism and antibody fluids
- antistatic treatment
- great air permeability
- high strength and reliability
- reinforced pad 38 g (PP + PE) on the chest (50x70 cm) and sleeve (30 cm)
- set-in sleeve

Material: SMMS (spunbond + 2 meltblown + spunbond non wovens).

Weight per piece: 120 ± 2g

Sterile gowns are individually packed.

Sterilized by oxide Ethylene.

Multilingual instructions: GB, FR, IT, ES, PT, DE, GR, SE, RO, PL, Arabic.