

TAPED SEAM BOOT COVERS - Type 4B-5B-6B



Product Code	21092
Unit of sale	box of 20 pcs
Minimum order	1
Type	PPE
Class	III

UK-REP	Not required
CH-REP	No
EAN/UPC	8023279210927

Description

TAPED SEAM BOOT COVERS - TYPE 4B-5B-6B - box of 20

- Size: 45 x 36 x 26 cm
- TYPE 4B-5B-6B EN 14126

Microporous disposable boot covers for general purposes as: cleanrooms, hospitals, labs, food industry, medical and pharmaceutical.

Optimized body fit and lightweight (60 g) to improve comfort and safety.

The PP material provides high level of comfort and protection against solid particles such as dust.

Features:

- PP + PE, microporous
- breathable and dust proof
- liquid penetration resistance
- barrier to fine particulates
- air and water vapor permeable to reduce risk of heat stress and fiber contamination
- with ties, elastic open

Multilingual instructions: GB, FR, IT, ES, PT, DE, GR, SE.