

DISPOSABLE GOWNS 25g/m2 115x137 cm - size M - non sterile



Product Code	21131
Unit of sale	box of 100 pcs
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	No
RDM (NSIS)	2527750
CND	T0205
EAN/UPC	6973866559811
UMDNS	15037
PARAF	983363173

Description

WATERPROOF DISPOSABLE GOWNS PP 25 g

- NON STERILE
- Size: Medium 115 x 137 cm

It can be used in clinics, patient rooms, laboratory, check point, community...
Not suitable for ICU and operating room.

- collar with two ties
- cuffs with elastic
- long sleeve
- waist with ties

Packing size: 52 x 33 x 30 cm

Gross weight: approx. 6.8 kg

Multilingual instructions: GB, FR, IT, ES, PT, DE, GR, SE, RO, PL, Arabic.

Sold in a box of 100 (10 pouches of 10 pieces).