

SUCTION CATHETER ch/fr 14 - 90 cm - sterile



Product Code	22125
Unit of sale	box of 100 pcs
Minimum order	1
Type	Medical device
Class	II A

UK-REP	No
CH-REP	Yes

RDM (NSIS)	2834777
CND	A06010103
EAN/UPC	8023279221251
PARAF	977810112

Description

PVC DEHP FREE* STERILE SUCTION CATHETERS with conical adapter

*DEHP is an additive used in the plastic industry for PVC which could be harmful to human reproductive system. Latex is totally DEHP free

Suction catheter ch/fr 14 - dark green
Length: 90 cm

Special technology: inner Ø does not shrink at edge

Multilingual box and manual: GB, FR, IT, ES, PT, DE, NL, PL, CZ, SK, HU, EE, LT, LV, DK, SE, FI, RO, BG, GR, HR, SI