



|               |                   |
|---------------|-------------------|
| Product Code  | 21402             |
| Unit of sale  | 1 pc              |
| Minimum order | 1                 |
| Type          | No medical device |
| Class         |                   |

|         |               |
|---------|---------------|
| EAN/UPC | 8023279214024 |
| PARAF   | 977809591     |

## Description

### DOCTOR'S WHITE COATS 60% COTTON - 40% POLYESTER

- WOMAN - With 5 buttons
- Size: (Europe: 42-44) (USA: S) (France: 38-40)
- 60% cotton / 40% polyester

Professional and stylish line with all sizes for men and women, made of 60% high quality cotton and 40% polyester.  
Doctor's white coats have 2 side pockets and 1 breast pocket.  
Man and woman model have buttons, unisex model has hidden press studs and elasticated cuffs.  
Weight 200 g/m<sup>2</sup>

## Technical Specifications

Composition: 60% cotton - 40% polyester  
Weight: 200 g/m<sup>2</sup>  
Yarn count: warp: Ne 20 ± 2 / weft: Ne 18 ± 2  
Construction: warp: 40 ± 3 / weft: 19 ± 3  
Breaking strength - warp: 100 kg ±3%  
Breaking strength - weft: 39 kg ±3%  
Fastness to washing: 4  
Treatments and finishing: crease resist  
Dimensional stability - washing at 90°C: 3%