



Product Code	21416
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

EAN/UPC	8023279214161
PARAF	977809678

## Description

### DOCTOR'S WHITE COATS 60% COTTON - 40% POLYESTER

- MAN - With 4 buttons
- Size: (Europe: 58-60) (USA: XXL) (France: 54-56)
- 60% cotton / 40% polyester

Professional and stylish line with all sizes for men and women, made of 60% high quality cotton and 40% polyester.

Doctor's white coats have 2 side pockets and 1 breast pocket.

Man and woman model have buttons, unisex model has hidden press studs and elasticated cuffs.

Weight 200 g/m<sup>2</sup>

## Technical Specifications

Composition: 60% cotton - 40% polyester

Weight: 200 g/m<sup>2</sup>

Yarn count: warp: Ne 20 ± 2 / weft: Ne 18 ± 2

Construction: warp: 40 ± 3 / weft: 19 ± 3

Breaking strength - warp: 100 kg ±3%

Breaking strength - weft: 39 kg ±3%

Fastness to washing: 4

Treatments and finishing: crease resist

Dimensional stability - washing at 90°C: 3%