



Product Code	21536
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

EAN/UPC	8023279215366
---------	---------------

## Description

TROUSERS - UNISEX - 60% COTTON - 40% POLYESTER

- Size: XXL
- Colour: blue navy
- UNISEX
- Elasticated waist
- 2 front and 1 back pockets

Professional and stylish line with all sizes for men and women, made of 60% high quality cotton and 40% polyester. Stock in 5 colours, other 22 colours available in 60 days.

## Technical Specifications

Composition: 60% cotton - 40% polyester  
 Weight: 200 g/m2 (white), 185 g/m2 (coloured)  
 Yarn count: warp: Ne 20 ± 2 weft: Ne 18 ± 2  
 Construction: warp: 40 ± 3 weft: 19 ± 3  
 Breaking strength - warp: 100 kg ±3%  
 Breaking strength - weft: 39 kg ±3% (white) 49 kg ±3% (coloured)  
 Fastness to washing: 4  
 Treatments and finishing: crease resist  
 Dimensional stability - washing at 90°C: 3%