

RIESTER DERMATOSCOPE LENS for 32150 - 13270



Product Code	32154
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

EAN/UPC	4045396211988
---------	---------------

Description

RCS-100 accessories: Dermatoscope lens with digital ruler for measurements

Riester code: 13270

Technical Specifications

Focal length/width: 2.2
Polarization: dermis / epidermis
Field of view: 30mm Ø
Object Distance: 0 mm
Magnification optical: 2x / digital: 4x / dispalyed: 2x-8x
Photography source: natural light LED
LED color temperature: 4000k
Size - weight: 62.1x62.1x36 mm - 108.5 g