

MIC Wi-Fi & USB POLARIZED+UV+WHITE LED DERMATOSCOPE with software



Product Code	32177
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	No

RDM (NSIS)	2184077
CND	Z12040118
EAN/UPC	8050713880132

Description

WI-FI - USB DIGITAL POLARIZED+UV+WHITE LED DERMATOSCOPE with software

Dermatoscope with polarized filter is a very innovative handheld colour digital device with USB and WI-FI transmission. Suitable for skin and scalp. Small, easy to operate and portable.
3 types of lighting: white, UV and white polarized LEDs.

Magnification:

60x and 200x with 7 mm spacer ring or closed cap
30x with transparent tube

MDR 2017/745/UE - Medical Device Class I

Mic-Fi Derma software for Windows (GB, IT) and MAC (GB) to capture photos, compare images, history and reporting, data base for patients and selection of the point where analysis is carried out.

Technical Specifications

Operating mode: continuous
Frame rate: 15 frame/sec.
Camera: Colour CMOS 1/4"
Image format: bmp, jpg
LEDs: 4 white, 4 UV and 4 polarized adjustable LEDs
USB: USB 2.0, 2 mt
Battery: Li-ion
Working time: 2 hours
Lens: high definition microscopy lens
Power supply: DC 5.0 V/1A
Consumption: Max 2.5 W
Resolution: 1.3 Mpixel
Size - weight: 36x142 mm - 88 g

Standard accessories

USB cable
Ruler
1 Bag L + 1 bag S
USB AC Adapter
Contact closed cap
Quick guide + manual (GB, FR, IT, ES, DE)
USB foot pedal

