



| | |
|---------------|----------------|
| Product Code | 21851 |
| Unit of sale | box of 100 pcs |
| Minimum order | 1 |
| Type | PPE |
| Class | III |

| | |
|---------|---------------|
| UK-REP | Not required |
| EAN/UPC | 5906615011798 |
| PARAF | 987476443 |

Description

SANTEX® POWDERED - ambidextrous

Small size - non sterile
Creamy colour

Non-sterile, powdered, examination and protective glove for single use.
Made of natural rubber latex.
Smooth external surface guarantees a tactile sensitivity during the work.
They minimize hand sweating.
Enhanced tactile sensitivity.
AQL=1.0

Multilingual box and instructions: GB, FR, IT, ES, PT, DE, PL, GR, SE, NO, FI, DK, SK, CZ, BG, RO, HU, HR, RS, SL, KZ, LT, LV, RU, TR.

Technical Specifications

External surface: Smooth
Cuff: beaded
MD (Medical Device): Class I
PPE (Personal protective equipment): Category III