

FOLEY 2-WAY LATEX CATHETER ch/fr 20 - ball 5-10 ml - sterile



Product Code	22058
Unit of sale	box of 10 pcs
Minimum order	1
Type	Medical device
Class	II A

UK-REP	No
CH-REP	Yes
RDM (NSIS)	2842241
CND	U010201
EAN/UPC	8023279220582
PARAF	977809805

Description

LATEX STERILE 100% SILICONE COATED FOLEY CATHETERS - SINGLE USE

- ch/fr 20

Two-way catheters - balloon 5-10 ml - length 40 cm - hard valve
Ball 5-10 ml means that the ball could be from 5 ml to 10 ml

Multilingual box and manual: GB, FR, IT, ES, PT, DE, NL, PL, CZ, SK, HU, EE, LT, LV, DK, SE, FI, RO, BG, GR, HR, SI