

# FOLEY 2-WAY LATEX CATHETER ch/fr 22 - ball 30 ml - sterile



Product Code	22069
Unit of sale	box of 50 pcs
Minimum order	1
Type	Medical device
Class	II A

UK-REP	No
CH-REP	Yes
RDM (NSIS)	2842241
CND	U010201
EAN/UPC	8023279220698
PARAF	977809870

## Description

LATEX STERILE 100% SILICONE COATED FOLEY CATHETERS - SINGLE USE

- ch/fr 22

Two-way catheters - balloon 30 ml - length 40 cm - hard valve

Multilingual box and manual: GB, FR, IT, ES, PT, DE, NL, PL, CZ, SK, HU, EE, LT, LV, DK, SE, FI, RO, BG, GR, HR, SI