

CATHETER LUBRICANT GEL 12 ml - sterile



Product Code	22104
Unit of sale	box of 25 pcs
Minimum order	1
Type	Medical device
Class	III

UK-REP	No
CH-REP	No

RDM (NSIS)	2971718
CND	V9099
EAN/UPC	8698712453737
PARAF	977809983
GMDN	60796

Description

LIDO C CATHETER LUBRICANT GEL 12 ml - with Lidocaine - sterile - box of 25

Sterile catheter gels are sterile, water-soluble, non-irritating lubricating gels with supportive local anesthetic effect with 2% Lidocaine and antiseptic properties with 0.05% Chlorhexidine. They are used to facilitate insertion of catheters and other medical instruments (catheterization, cystoscopy) into the urethra and bladder. They reduce the risk of urethral trauma during catheter insertion and possibility of infection by reducing friction. Latex free and hypoallergenic. Effective and long-lasting lubrication. It can be applied easily. It has ideal viscosity. Class III

After Gamma sterilization it may change the colour over time from transparent to yellowish or pinkish.