

SUCTION CATHETER ch/fr 6 - 50 cm - sterile



| | |
|---------------|----------------|
| Product Code | 22110 |
| Unit of sale | box of 100 pcs |
| Minimum order | 1 |
| Type | Medical device |
| Class | II A |

| | |
|--------|-----|
| UK-REP | No |
| CH-REP | Yes |

| | |
|------------|---------------|
| RDM (NSIS) | 2834777 |
| CND | A06010103 |
| EAN/UPC | 8023279221107 |
| PARAF | 977809995 |

Description

PVC DEHP FREE* STERILE SUCTION CATHETERS with conical adapter

*DEHP is an additive used in the plastic industry for PVC which could be harmful to human reproductive system. Latex is totally DEHP free

Suction catheter ch/fr 6 - light green
Length: 50 cm

Special technology: inner Ø does not shrink at edge

Multilingual box and manual: GB, FR, IT, ES, PT, DE, NL, PL, CZ, SK, HU, EE, LT, LV, DK, SE, FI, RO, BG, GR, HR, SI