

# SUCTION CATHETER ch/fr 8 - 50 cm - sterile



Product Code	22112
Unit of sale	box of 100 pcs
Minimum order	1
Type	Medical device
Class	II A

UK-REP	No
CH-REP	Yes

RDM (NSIS)	2834777
CND	A06010103
EAN/UPC	8023279221121
PARAF	977810009

## Description

PVC DEHP FREE\* STERILE SUCTION CATHETERS with conical adapter

\*DEHP is an additive used in the plastic industry for PVC which could be harmful to human reproductive system. Latex is totally DEHP free

Suction catheter ch/fr 8 - light blue  
Length: 50 cm

Special technology: inner Ø does not shrink at edge

Multilingual box and manual: GB, FR, IT, ES, PT, DE, NL, PL, CZ, SK, HU, EE, LT, LV, DK, SE, FI, RO, BG, GR, HR, SI