

# ETHICON ETHILON MONOFILAMENT SUTURES - gauge 4/0 needle 13 mm



Product Code	22294
Unit of sale	box of 36 pcs
Minimum order	1
Type	Medical device
Class	III

UK-REP	No
CH-REP	No

RDM (NSIS)	2561759
CND	H010201010101
GTIN	10705031008632

## Description

ETHILON - non absorbable sterile monofilament sutures

- Ethicon code: 699H
- Gauge: 4/0
- Metric: 1.5
- Colour: black
- Circle: 3/8
- Needle length: 13 mm
- Needle: P-3
- Suture length: 45 cm
- Minimum order: box of 36 pieces

Ethilon is indicated for use in general soft tissue approximation and/or ligation, including cardiovascular, ophthalmic and neurosurgical procedures.