

ETHICON ETHILON MONOFILAMENT SUTURES - gauge 4/0 needle 16 mm



Product Code	22298
Unit of sale	box of 36 pcs
Minimum order	1
Type	Medical device
Class	III

UK-REP	No
CH-REP	No

RDM (NSIS)	2559897
CND	H010201010101
EAN/UPC	8023279222982
GMDN	38000

Description

ETHILON - non absorbable sterile monofilament sutures

- Ethicon code: 1162H
- Gauge: 4/0
- Metric: 1.5
- Colour: black
- Circle: 3/8
- Needle length: 16 mm
- Needle: FS-3
- Suture length: 45 cm
- Minimum order: box of 36 pieces

Non absorbable surgical sutures are indicated for use in general soft tissue approximation and/or ligation included use in cardiovascular, ophthalmic and neurosurgical procedures.