

# ETHICON ETHILON MONOFILAMENT SUTURES - gauge 3/0 needle 19 mm



Product Code	22304
Unit of sale	box of 36 pcs
Minimum order	1
Type	Medical device
Class	III
UK-REP	No
CH-REP	No
RDM (NSIS)	2561751
CND	H010201010101
EAN/UPC	8023279223040
GMDN	38000

## Description

ETHILON - non absorbable sterile monofilament sutures

- Ethicon code: EH7506H
- Gauge: 3/0
- Metric: 2.0
- Colour: black
- Circle: 3/8
- Needle length: 19 mm
- Needle: FS-2
- Suture length: 45 cm
- Minimum order: box of 36 pieces

Non absorbable surgical sutures are indicated for use in general soft tissue approximation and/or ligation included use in cardiovascular, ophthalmic and neurosurgical procedures.