

ETHICON PERMA-HAND SILK SUTURES - gauge 3/0 needle 17 mm - braided



Product Code	22317
Unit of sale	box of 36 pcs
Minimum order	1
Type	Medical device
Class	III

UK-REP	No
CH-REP	No

RDM (NSIS)	2471108
CND	H010202020301
GTIN	30705031100682
GMDN	13910

Description

PERMA-HAND - non absorbable sterile silk sutures, braided

- Ethicon code: K872
- Gauge: 3/0
- Metric: 2.0
- Colour: black
- Circle: 1/2
- Needle length: 17 mm
- Needle: RB-1
- Suture length: 75 cm
- Minimum order: box of 36 pieces

Non absorbable surgical sutures are indicated for use in general soft tissue approximation and/or ligation included use in cardiovascular, ophthalmic and neurosurgical procedures.