

ETHICON PERMA-HAND SILK SUTURES - gauge 3/0 needle 24 mm - braided



Product Code	22319
Unit of sale	box of 36 pcs
Minimum order	1
Type	Medical device
Class	III

UK-REP	No
CH-REP	No

RDM (NSIS)	2471030
CND	H010202020301
GTIN	30705031090051
GMDN	13910

Description

PERMA-HAND - non absorbable sterile silk sutures, braided

- Ethicon code: 684H
- Gauge: 3/0
- Metric: 2.0
- Colour: black
- Circle: 3/8
- Needle length: 24 mm
- Needle: FS-1
- Suture length: 45 cm
- Minimum order: box of 36 pieces

Non absorbable surgical sutures are indicated for use in general soft tissue approximation and/or ligation included use in cardiovascular, ophthalmic and neurosurgical procedures.