

# ETHICON PROLENE BLUE MONOFILAMENT SUTURES - gauge 4/0 needle 19 mm FS-2 suture 75 cm



Product Code	22320
Unit of sale	box of 36 pcs
Minimum order	1
Type	Medical device
Class	III
UK-REP	No
CH-REP	No
RDM (NSIS)	2558542
CND	H010201010401
EAN/UPC	8023279223200
GMDN	13909

## Description

PROLENE - blue pigmented monofilament sutures, non absorbable sterile

- Ethicon code: EH7150H
- Gauge: 4/0
- Metric: 1.5
- Colour: blue
- Circle: 3/8
- Needle length: 19 mm
- Needle: FS-2
- Suture length: 75 cm
- Minimum order: box of 36 pieces

Non absorbable surgical sutures are indicated for use in general soft tissue approximation and/or ligation included use in cardiovascular, ophthalmic and neurosurgical procedures.