

ETHICON PROLENE BLUE MONOFILAMENT SUTURES

- gauge 3/0 needle 19 mm FS-3



Product Code	22324
Unit of sale	box of 36 pcs
Minimum order	1
Type	Medical device
Class	III

UK-REP	No
CH-REP	No
RDM (NSIS)	2558678
CND	H010201010401
EAN/UPC	8023279223248
GMDN	13909

Description

PROLENE - blue pigmented monofilament sutures, non absorbable sterile

- Ethicon code: 8665H
- Gauge: 3/0
- Metric: 2.0
- Colour: blue
- Circle: 3/8
- Needle lenght: 19 mm
- Needle: FS-3
- Suture lenght: 45 cm
- Minimum order: box of 36 pieces

Non absorbable surgical sutures are indicated for use in general soft tissue approximation and/or ligation included use in cardiovascular, ophthalmic and neurosurgical procedures.