

# ETHICON PROLENE BLUE MONOFILAMENT SUTURES

## - gauge 4/0 needle 19 mm PS-2



Product Code	22325
Unit of sale	box of 36 pcs
Minimum order	1
Type	Medical device
Class	III

UK-REP	No
CH-REP	No
RDM (NSIS)	1425909
CND	H010201010401
EAN/UPC	8023279223255
GMDN	13909

### Description

PROLENE - blue pigmented monofilament sutures, non absorbable sterile

- Ethicon code: 8682H
- Gauge: 4/0
- Metric: 1.5
- Colour: blue
- Circle: 3/8
- Needle length: 19 mm
- Needle: PS-2
- Suture length: 45 cm
- Minimum order: box of 36 pieces

Non absorbable surgical sutures are indicated for use in general soft tissue approximation and/or ligation included use in cardiovascular, ophthalmic and neurosurgical procedures.