

ETHICON PROLENE BLUE MONOFILAMENT SUTURES

- gauge 4/0 needle 13 mm P-3



Product Code	22329
Unit of sale	box of 36 pcs
Minimum order	1
Type	Medical device
Class	III

UK-REP	No
CH-REP	No

RDM (NSIS)	2558662
CND	H010201010401
EAN/UPC	8023279223293
GMDN	13909

Description

PROLENE - blue pigmented monofilament sutures, non absorbable sterile

- Ethicon code: 8699H
- Gauge: 4/0
- Metric: 1.5
- Colour: blue
- Circle: 3/8
- Needle length: 13 mm
- Needle: P-3
- Suture length: 45 cm
- Minimum order: box of 36 pieces

Non absorbable surgical sutures are indicated for use in general soft tissue approximation and/or ligation included use in cardiovascular, ophthalmic and neurosurgical procedures.