

ETHICON PROLENE BLUE MONOFILAMENT SUTURES

- gauge 6/0 needle 13 mm P-3 Prime MP



Product Code	22330
Unit of sale	box of 36 pcs
Minimum order	1
Type	Medical device
Class	III

UK-REP	No
CH-REP	No

RDM (NSIS)	2558665
CND	H010201010401

Description

PROLENE - blue pigmented monofilament sutures, non absorbable sterile

- Ethicon code: 8695H
- Gauge: 6/0
- Metric: 0.7
- Colour: blue
- Circle: 3/8
- Needle length: 13 mm
- Needle: P-3 Prime MP
- Suture length: 45 cm
- Minimum order: box of 36 pieces

Prolene is indicated for vessels and soft tissues, sutures are particularly suitable for cardiovascular, neurosurgical and ophthalmic procedures.