

ETHICON MONOCRYL ABSORBABLE SUTURES - gauge 4/0 needle 16 mm



Product Code	22351
Unit of sale	box of 12 pcs
Minimum order	1
Type	Medical device
Class	III
UK-REP	No
CH-REP	No
RDM (NSIS)	2541978
CND	H010101010301
EAN/UPC	8023279223514
GMDN	45814

Description

MONOCRYL - absorbable sterile surgical sutures

Complete absorption in 91-119 days

- Ethicon code: MCP3205G
- Gauge: 4/0
- Metric: 1.5
- Colour: undyed
- Circle: 3/8
- Needle lenght: 16 mm
- Needle: PS-3
- Suture lenght: 70 cm
- Minimum order: box of 12 pieces

Monocryl is indicated for use in general soft tissue approximation and/or ligation, but not for cardio, neuro, micro and ophthalmic surgery.