

ETHICON VICRYL PLUS ABSORBABLE SUTURES - gauge 3/0 needle 22 mm - braided



Product Code	22364
Unit of sale	box of 36 pcs
Minimum order	1
Type	Medical device
Class	III

UK-REP	No
CH-REP	No

RDM (NSIS)	2487064
CND	H010101040201
EAN/UPC	8023279223644
GMDN	45401

Description

VICRYL PLUS BRAIDED - absorbable sterile surgical sutures

Complete absorption in 56-70 days

- Ethicon code: VCP311H
- Gauge: 3/0
- Metric: 2.0
- Colour: violet
- Circle: 3/8
- Needle length: 22 mm
- Needle: SH-1
- Suture length: 70 cm
- Minimum order: box of 36 pieces

vicryl plus is indicated for use in general tissue approximation and/or ligation except for ophthalmic, neuro, cardiovascular.