

ETHICON VICRYL RAPID ABSORBABLE SUTURES - gauge 3/0 needle 17 mm - braided



| | |
|---------------|----------------|
| Product Code | 22370 |
| Unit of sale | box of 36 pcs |
| Minimum order | 1 |
| Type | Medical device |
| Class | III |

| | |
|------------|----------------|
| UK-REP | No |
| CH-REP | No |
| RDM (NSIS) | 2436988 |
| CND | H0101010299 |
| EAN/UPC | 8023279223705 |
| GTIN | 30705031109630 |
| GMDN | 17471 |

Description

VICRYL RAPID BRAIDED - absorbable sterile surgical sutures

Complete absorption in 42 days

- Ethicon code: V76H
- Gauge: 3/0
- Metric: 2.0
- Colour: undyed
- Needle: J-1
- Suture length: 45 cm
- Minimum order: box of 36 pieces

Vicryl Rapid Braided indicated only for: use in superficial soft tissue approximation of the skin and mucosa, where only short term (7-10 days) wound support is required.