

# ETHICON VICRYL RAPID ABSORBABLE SUTURES - gauge 5/0 needle 19 mm - braided



Product Code	22376
Unit of sale	box of 36 pcs
Minimum order	1
Type	Medical device
Class	III

RDM (NSIS)	2436963
CND	H0101010299
GTIN	30705031106523
GMDN	17471

## Description

VICRYL RAPID BRAIDED - absorbable sterile surgical sutures

Complete absorption in 42 days

- Ethicon code: V2910H
- Gauge: 5/0
- Metric: 1.0
- Colour: undyed
- Circle: 3/8C
- Needle length: 19 mm
- Needle: FS-2
- Suture length: 45 cm
- Minimum order: box of 36 pieces

Vicryl Rapid Braided indicated only for: use in superficial soft tissue approximation of the skin and mucosa, where only short term (7-10 days) wound support is required.