

# SWANN-MORTON CARBON STEEL STITCH CUTTER BLADES - sterile



Product Code	22420
Unit of sale	box of 100 pcs
Minimum order	1
Type	Medical device
Class	I S

UK-REP	Yes
CH-REP	Yes

RDM (NSIS)	2370294
CND	V0102
EAN/UPC	8023279224207
GMDN	16224

## Description

### CARBON STEEL STITCH CUTTER BLADE - STERILE

- Carbon steel stitch cutter blade

Multilingual box: GB, FR, IT, ES, DE, DK, NL, GR, PT, PL, FI, RU, NO, SE.

Swann-Morton blades and scalpels are renowned their exceptional quality, precision, consistency and reliability. Swann-Morton manufactures in the UK using high quality steel that gives exceptional cutting performance. They operate stringent quality control procedures which ensure unrivalled blade consistency and reliability. This is why Swann-Morton is the preferred choice of many medical professionals worldwide.

Made in the UK.