

SWANN-MORTON STAINLESS STEEL HANDLE N. 4 without CE



Product Code	22424
--------------	-------

Unit of sale	box of 10 pcs
--------------	------------------

Minimum order	1
---------------	---

Type	No medical device
------	-------------------

Class	
-------	--

UK-REP	Not required
--------	--------------

EAN/UPC	5033955010348
---------	---------------

Description

SWANN-MORTON S/S HANDLE

- Graduated S/S n° 4 handle
- Compatibility with blade No.: 20, 21, 22, 23, 24
- Without CE
- Box of 10

Multilingual box: GB, FR, IT, ES, DE, DK, NL, GR, PT, PL, FI, RU, NO, SE.

Swann-Morton blades and scalpels are renowned their exceptional quality, precision, consistency and reliability. Swann-Morton manufactures in the UK using high quality steel that gives exceptional cutting performance. They operate stringent quality control procedures which ensure unrivalled blade consistency and reliability. This is why Swann-Morton is the preferred choice of many medical professionals worldwide.

Made in the UK.