

# SWANN-MORTON STAINLESS STEEL HANDLE N. 3 with CE and UKCA



Product Code	22425
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I R

UK-REP	Yes
CH-REP	Yes

RDM (NSIS)	2320379
CND	L010103
EAN/UPC	8023279224252
GMDN	12235

## Description

### SWANN-MORTON S/S HANDLE

- Graduated S/S n° 3 handle
- Compatibility with blade No.: 10, 11, 12, 15
- With CE

Multilingual box: GB, FR, IT, ES, DE, DK, NL, GR, PT, PL, FI, RU, NO, SE.

Swann-Morton blades and scalpels are renowned their exceptional quality, precision, consistency and reliability. Swann-Morton manufactures in the UK using high quality steel that gives exceptional cutting performance. They operate stringent quality control procedures which ensure unrivalled blade consistency and reliability. This is why Swann-Morton is the preferred choice of many medical professionals worldwide.

Made in the UK.