


SWANN-MORTON STAINLESS STEEL HANDLE N. 4 with CE and UKCA



Product Code	22426
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I R

UK-REP	Yes
CH-REP	Yes
RDM (NSIS)	2320371
CND	L010103
EAN/UPC	8023279224269
GMDN	12235

Description

SWANN-MORTON S/S HANDLE

- Graduated S/S n° 4 handle
- Compatibility with blade No.: 20, 21, 22, 23, 24
- With CE

Multilingual box: GB, FR, IT, ES, DE, DK, NL, GR, PT, PL, FI, RU, NO, SE.

Swann-Morton blades and scalpels are renowned their exceptional quality, precision, consistency and reliability. Swann-Morton manufactures in the UK using high quality steel that gives exceptional cutting performance. They operate stringent quality control procedures which ensure unrivalled blade consistency and reliability. This is why Swann-Morton is the preferred choice of many medical professionals worldwide.

Made in the UK.