

# SWANN-MORTON SCALPELS WITH STAINLESS STEEL BLADE N. 23 - sterile



Product Code	22438
Unit of sale	box of 10 pcs
Minimum order	1
Type	Medical device
Class	II A

UK-REP	Yes
CH-REP	Yes

RDM (NSIS)	2370305
CND	V010102
EAN/UPC	5033955005108
GTIN	05033955005108
GMDN	47569

## Description

### SWANN-MORTON SCALPELS WITH S/S BLADE

- Disposable scalpels - blade 23 - sterile

Plastic handle scalpels with stainless steel blades, individually packed in sterilized envelope in box of 10 pcs.

Multilingual box: GB, FR, IT, ES, DE, DK, NL, GR, PT, PL, FI, RU, NO, SE.

Swann-Morton blades and scalpels are renowned their exceptional quality, precision, consistency and reliability. Swann-Morton manufactures in the UK using high quality steel that gives exceptional cutting performance. They operate stringent quality control procedures which ensure unrivalled blade consistency and reliability. This is why Swann-Morton is the preferred choice of many medical professionals worldwide.

Made in the UK.