

# SWANN-MORTON RETRACTABLE SCALPELS WITH S/S BLADE N. 23 - sterile



Product Code	22448
Unit of sale	box of 25 pcs
Minimum order	1
Type	Medical device
Class	II A

UK-REP	Yes
CH-REP	Yes
RDM (NSIS)	2543031
CND	V010101
EAN/UPC	8023279224481
GMDN	47569

## Description

### SWANN-MORTON SAFETY SCALPELS WITH RETRACTABLE STAINLESS STEEL BLADE

- Disposable scalpels - blade 23 - sterile

Multilingual box: GB, FR, IT, ES, DE, DK, NL, GR, PT, PL, FI, RU, NO, SE.

Swann-Morton blades and scalpels are renowned their exceptional quality, precision, consistency and reliability. Swann-Morton manufactures in the UK using high quality steel that gives exceptional cutting performance. They operate stringent quality control procedures which ensure unrivalled blade consistency and reliability. This is why Swann-Morton is the preferred choice of many medical professionals worldwide.

Made in the UK.