

SWANN-MORTON STITCH CUTTER WITH RETRACTABLE S/S BLADE - sterile



Product Code	22450
Unit of sale	box of 25 pcs
Minimum order	1
Type	Medical device
Class	I S

UK-REP	Yes
CH-REP	Yes

RDM (NSIS)	1420670
CND	V0101
EAN/UPC	8023279224504
GMDN	16224

Description

STITCH CUTTER WITH RETRACTABLE S/S BLADE

- Stitch cutter with retractable blade - sterile

Multilingual box: GB, FR, IT, ES, DE, DK, NL, GR, PT, PL, FI, RU, NO, SE.

Swann-Morton blades and scalpels are renowned their exceptional quality, precision, consistency and reliability. Swann-Morton manufactures in the UK using high quality steel that gives exceptional cutting performance. They operate stringent quality control procedures which ensure unrivalled blade consistency and reliability. This is why Swann-Morton is the preferred choice of many medical professionals worldwide.

Made in the UK.