

DIAMANTINE DISPOSABLE SCALPELS WITH S/S BLADE N. 24 - sterile



Product Code	22468
Unit of sale	box of 10 pcs
Minimum order	1
Type	Medical device
Class	II A

UK-REP	Yes
CH-REP	No
RDM (NSIS)	2391016
CND	V010301
EAN/UPC	8023279224689

Description

DIAMANTINE™ EDGE HIGH PERFORMANCE SCALPELS - sterile

- Scalpels - blade 24
- Sterile

High-Performance Edge blade for extreme and durable sharpness.

- sensory grip handle enhances cutting action and control of the line of incision
- each handle has a graduated 6 cm ruler to help determine incision length (indicative only)
- protective sheath covers blade to help prevent accidental cuts
- latex and phthalates free

Sensory grip handle.
Stainless steel Diamantine Edge high performance coated blade.
Hardness over 2,500 Vickers

Size is clearly marked on blade and handle for easy identification.